



INNOVATIVE
ORTHO MEDICAL

Product
Catalogue



INNOVATIVE
ORTHO MEDICAL

Trauma Plate
System

Middle S-clavicle Plate
 (Use 3.5mm Cortical Screw, 4.0mm
 Cancellous Screw)



8 holes L	S4100-0101
8 holes R	S4100-0102

Distal S-clavicle Plate
 (Use 3.5mm Cortical Screw, 4.0mm
 Cancellous Screw)



4 holes L	S4100-0201
6 holes L	S4100-0202
8 holes L	S4100-0203
10 holes L	S4100-0204
4 holes R	S4100-0205
6 holes R	S4100-0206
8 holes R	S4100-0207
10 holes R	S4100-0208



Distal S-clavicle Plate-II
(Use 3.5mm Cortical Screw, 4.0mm
Cancellous Screw)



7 holes L	S4100-0301
8 holes L	S4100-0302
9 holes L	S4100-0303
10 holes L	S4100-0304
7 holes R	S4100-0305
8 holes R	S4100-0306
9 holes R	S4100-0307
10 holes R	S4100-0308

Clavicle Claw Plate
(Use 3.5mm Cortical Screw, 4.0mm
Cancellous Screw)



4 holes L	S4100-0401
5 holes L	S4100-0402
6 holes L	S4100-0403
4 holes R	S4100-0404
5 holes R	S4100-0405
6 holes R	S4100-0406



LC-DCP Unla and Radius Plate
(Use 3.5mm Cortical Screw)



5 holes	S4100-0601
6 holes	S4100-0602
7 holes	S4100-0603
8 holes	S4100-0604
9 holes	S4100-0605
10 holes	S4100-0606

One-third Tubular Plate
(Use 3.5mm Cortical Screw)



4 holes	S4100-0701
5 holes	S4100-0702
6 holes	S4100-0703
7 holes	S4100-0704
8 holes	S4100-0705

DCP Humeral Plate
(Use 3.5mm Cortical Screw)




4 holes	S4100-0801
5 holes	S4100-0802
6 holes	S4100-0803
7 holes	S4100-0804
8 holes	S4100-0805



9 holes	S4100-0806
10 holes	S4100-0807
11 holes	S4100-0808
12 holes	S4100-0809



LC-DCP Humeral Plate
(Use 3.5mm Cortical Screw)




4 holes	S4100-0901
5 holes	S4100-0902
6 holes	S4100-0903
7 holes	S4100-0904
8 holes	S4100-0905
9 holes	S4100-0906
10 holes	S4100-0907
11 holes	S4100-0908
12 holes	S4100-0909

LC-DCP Humerus Plate-II
(Use 4.5mm Cortical Screw)



5 holes	S4100-1001
6 holes	S4100-1002
7 holes	S4100-1003
8 holes	S4100-1004
9 holes	S4100-1005
10 holes	S4100-1006



Small T-type Curved Flate Plate
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)




3 holes L	S4100-1101
4 holes L	S4100-1102
5 holes L	S4100-1103
3 holes R	S4100-1104
4 holes R	S4100-1105
5 holes R	S4100-1106

Distal Radius Medial Plate
(Metacarpus Medial Type)
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)



3 holes	S4100-1201
4 holes	S4100-1202
5 holes	S4100-1203
6 holes	S4100-1204

Distal Radius Back Plate
(Use 2.5mm Cortical Screw)



L	S4100-1401
R	S4100-1402



Distal Radius Medial Plate
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)



3 holes L	S4100-1301
4 holes L	S4100-1302
5 holes L	S4100-1303
6 holes L	S4100-1304
3 holes R	S4100-1305
4 holes R	S4100-1306
5 holes R	S4100-1307
6 holes R	S4100-1308

Distal Radius/Fibula Plate
(Use 3.5mm Cortical Screw,3.5mm
Cancellous Screw)



3 holes	S4100-1501
4 holes	S4100-1502
5 holes	S4100-1503
6 holes	S4100-1504



Olecranon Reconstruction Plate
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)



3 holes L	S4100-1801
4 holes L	S4100-1802
5 holes L	S4100-1803
6 holes L	S4100-1804
3 holes R	S4100-1805
4 holes R	S4100-1806
5 holes R	S4100-1807
6 holes R	S4100-1808

Humerus Condylus Plate
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)



3 holes	S4100-2001
4 holes	S4100-2002
5 holes	S4100-2003
6 holes	S4100-2004
7 holes	S4100-2005
8 holes	S4100-2006
9 holes	S4100-2007
10 holes	S4100-2008



Small T-type Curved Plate
(Use 3.5mm Cortical Screw, 3.5mm
Cancellous Screw)



3 holes

S4100-1601

4 holes

S4100-1602

5 holes

S4100-1603

Distal Radius Lateral Plate
(Use 3.5mm Cortical Screw, 3.5mm
Cancellous Screw)



3 holes

S4100-1604

4 holes

S4100-1605

5 holes

S4100-1606

6 holes

S4100-1607

Olecranon Plate
(Use 3.5mm Cortical Screw, 4.0mm
Cancellous Screw)



3 holes L

S4100-1701

4 holes L

S4100-1702

3 holes R

S4100-1703

4 holes R

S4100-1704



Olecranon Plate
 (Use 3.5mm Cortical Screw,4.0mm
 Cancellous Screw)



3 holes	S4100-1901
4 holes	S4100-1902
5 holes	S4100-1903
6 holes	S4100-1904

Humerus Condylus Plate
 (Use 3.5mm Cortical Screw,4.0mm
 Cancellous Screw)



3 holes	S4100-2001
4 holes	S4100-2002
5 holes	S4100-2003
6 holes	S4100-2004
7 holes	S4100-2005
8 holes	S4100-2006
9 holes	S4100-2007
10 holes	S4100-2008

Three-Flanged Plate
 (Use 4.0mm Cortical Screw,4.0mm
 Cancellous Screw)



3 holes	S4100-2301
4 holes	S4100-2302
5 holes	S4100-2303
6 holes	S4100-2304
7 holes	S4100-2305
8 holes	S4100-2306



Rib Plate



11*30	/	/	T4100-2401
11*45	/	/	T4100-2402
11*55	/	/	T4100-2403
15*30	/	/	T4100-2404
15*45	/	/	T4100-2405
15*55	/	/	T4100-2406
17*30	/	/	T4100-2407
17*45	/	/	T4100-2408
17*55	/	/	T4100-2409
19*30	/	/	T4100-2410
19*45	/	/	T4100-2411
19*55	/	/	T4100-2412

Y-olecranon Fossa Plate
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)



3 holes L	S4100-2501
4 holes L	S4100-2502
5 holes L	S4100-2503
6 holes L	S4100-2504
7 holes L	S4100-2505
8 holes L	S4100-2506
3 holes R	S4100-2507
4 holes R	S4100-2508
5 holes R	S4100-2509
6 holes R	S4100-2510
7 holes R	S4100-2511
8 holes R	S4100-2512



Y-shaped Humerus Condylus Plate
(Use 3.5mm Cortical Screw, 4.0mm
Cancellous Screw)



3 holes L	S4100-2601
4 holes L	S4100-2602
5 holes L	S4100-2603
6 holes L	S4100-2604
3 holes R	S4100-2605
4 holes R	S4100-2606
5 holes R	S4100-2607
6 holes R	S4100-2608

Distal Humeral ST Plate
(Use 3.5mm Cortical Screw, 4.0mm
Cancellous Screw)



3 holes L	S4100-2701
4 holes L	S4100-2702
5 holes L	S4100-2703
6 holes L	S4100-2704
3 holes R	S4100-2705
4 holes R	S4100-2706
5 holes R	S4100-2707
6 holes R	S4100-2708




Distal Humeral Lateral Plate
(Olecranon Type)
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)



12 holes L	S4100-2801
14 holes L	S4100-2802
16 holes L	S4100-2803
18 holes L	S4100-2804
12 holes R	S4100-2805
14 holes R	S4100-2806
16 holes R	S4100-2807
18 holes R	S4100-2808

Distal Radius Back Plate
(Use 2.5mm Cortical Screw)



L	S4100-1401
R	S4100-1402



Distal Humeral Medial Plate
(Reconstruction Type)
(Use 3.5mm Cortical Screw, 4.0mm
Cancellous Screw)



10 holes	S4100-2901
12 holes	S4100-2902
14 holes	S4100-2903

LC-DCP Tibial Plate
(Use 4.5mm Cortical Screw)



4 holes	S4100-3001
5 holes	S4100-3002
6 holes	S4100-3003
7 holes	S4100-3004
8 holes	S4100-3005
9 holes	S4100-3006
10 holes	S4100-3007
11 holes	S4100-3008
12 holes	S4100-3009



DCP Tibial Plate Narrow
(Use 4.5mm Cortical Screw)



4 holes	S4100-3101
5 holes	S4100-3102
6 holes	S4100-3103
7 holes	S4100-3104
8 holes	S4100-3105
9 holes	S4100-3106
10 holes	S4100-3107
11 holes	S4100-3108
12 holes	S4100-3109

T-shaped Plate
(Use 4.5mm Cortical Screw, 6.5mm
Cancellous Screw)



4 holes	S4100-3301
5 holes	S4100-3302
6 holes	S4100-3303
7 holes	S4100-3304
8 holes	S4100-3305



Reconstruction T-shaped Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



4 holes	S4100-3401
5 holes	S4100-3402
6 holes	S4100-3403
7 holes	S4100-3404
8 holes	S4100-3405

Distal Fibula Plate-I
(Use 3.5mm Cortical Screw,3.5mm
Cancellous Screw)



3 holes L	S4100-3501
4 holes L	S4100-3502
5 holes L	S4100-3503
6 holes L	S4100-3504
7 holes L	S4100-3505
8 holes L	S4100-3506
3 holes R	S4100-3507
4 holes R	S4100-3508
5 holes R	S4100-3509
6 holes R	S4100-3510
7 holes R	S4100-3511
8 holes R	S4100-3512



Distal Fibula Plate-II
(Use 3.5mm Cortical Screw,3.5mm
Cancellous Screw)




3 holes L	S4100-3601
4 holes L	S4100-3602
5 holes L	S4100-3603
6 holes L	S4100-3604
7 holes L	S4100-3605
8 holes L	S4100-3606
3 holes R	S4100-3607
4 holes R	S4100-3608
5 holes R	S4100-3609
6 holes R	S4100-3610
7 holes R	S4100-3611
8 holes R	S4100-3612

Distal Fibula Plate-III
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)



3 holes L	S4100-3701
4 holes L	S4100-3702
5 holes L	S4100-3703
6 holes L	S4100-3704
7 holes L	S4100-3705
8 holes L	S4100-3706
3 holes R	S4100-3707
4 holes R	S4100-3708
5 holes R	S4100-3709
6 holes R	S4100-3710
7 holes R	S4100-3711
8 holes R	S4100-3712



Distal Fibula Back Plate
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)




L-shaped Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



3 holes L	S4100-3801
4 holes L	S4100-3802
5 holes L	S4100-3803
6 holes L	S4100-3804
7 holes L	S4100-3805
8 holes L	S4100-3806
3 holes R	S4100-3807
4 holes R	S4100-3808
5 holes R	S4100-3809
6 holes R	S4100-3810
7 holes R	S4100-3811
8 holes R	S4100-3812

4 holes L	S4100-3901
5 holes L	S4100-3902
6 holes L	S4100-3903
7 holes L	S4100-3904
8 holes L	S4100-3905
4 holes R	S4100-3906
5 holes R	S4100-3907
6 holes R	S4100-3908
7 holes R	S4100-3909
8 holes R	S4100-3910



Reconstruction L-shaped Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



4 holes L	S4100-4001
5 holes L	S4100-4002
6 holes L	S4100-4003
7 holes L	S4100-4004
8 holes L	S4100-4005
4 holes R	S4100-4006
5 holes R	S4100-4007
6 holes R	S4100-4008
7 holes R	S4100-4009
8 holes R	S4100-4010



Tibial Platform Lateral Plate-I
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes L	S4100-4101
6 holes L	S4100-4102
7 holes L	S4100-4103
8 holes L	S4100-4104
9 holes L	S4100-4105
10 holes L	S4100-4106
11 holes L	S4100-4107
5 holes R	S4100-4108
6 holes R	S4100-4109
7 holes R	S4100-4110
8 holes R	S4100-4111
9 holes R	S4100-4112
10 holes R	S4100-4113
11 holes R	S4100-4114



Tibial Platform Lateral Plate-II
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes L	S4100-4201
6 holes L	S4100-4202
7 holes L	S4100-4203
8 holes L	S4100-4204
9 holes L	S4100-4205
10 holes L	S4100-4206
11 holes L	S4100-4207
5 holes R	S4100-4208
6 holes R	S4100-4209
7 holes R	S4100-4210
8 holes R	S4100-4211
9 holes R	S4100-4212
10 holes R	S4100-4213
11 holes R	S4100-4214

Tibial Platform Lateral Plate-III
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes L	S4100-4301
7 holes L	S4100-4302
9 holes L	S4100-4303
11 holes L	S4100-4304
5 holes R	S4100-4305
7 holes R	S4100-4306
9 holes R	S4100-4307
11 holes R	S4100-4308



Proximal Tibial Medial Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



4 holes L	S4100-4401
5 holes L	S4100-4402
6 holes L	S4100-4403
7 holes L	S4100-4404
8 holes L	S4100-4405
9 holes L	S4100-4406
10 holes L	S4100-4407
11 holes L	S4100-4408
4 holes R	S4100-4409
5 holes R	S4100-4410
6 holes R	S4100-4411
7 holes R	S4100-4412
8 holes R	S4100-4413
9 holes R	S4100-4414
10 holes R	S4100-4415
11 holes R	S4100-4416



Proximal Tibial Lateral Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



4 holes L	S4100-4501
5 holes L	S4100-4502
6 holes L	S4100-4503
7 holes L	S4100-4504
8 holes L	S4100-4505
9 holes L	S4100-4506
10 holes L	S4100-4507
11 holes L	S4100-4508
4 holes R	S4100-4509
5 holes R	S4100-4510
6 holes R	S4100-4511
7 holes R	S4100-4512
8 holes R	S4100-4513
9 holes R	S4100-4514
10 holes R	S4100-4515
11 holes R	S4100-4516



Distal Tibial Medial Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)




5 holes L	S4100-4601
6 holes L	S4100-4602
7 holes L	S4100-4603
8 holes L	S4100-4604
9 holes L	S4100-4605
10 holes L	S4100-4606
11 holes L	S4100-4607
5 holes R	S4100-4608
6 holes R	S4100-4609
7 holes R	S4100-4610
8 holes R	S4100-4611
9 holes R	S4100-4612
10 holes R	S4100-4613
11 holes R	S4100-4614



Distal Tibial Lateral Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes L	S4100-4701
6 holes L	S4100-4702
7 holes L	S4100-4703
8 holes L	S4100-4704
9 holes L	S4100-4705
10 holes L	S4100-4706
11 holes L	S4100-4707
5 holes R	S4100-4708
6 holes R	S4100-4709
7 holes R	S4100-4710
8 holes R	S4100-4711
9 holes R	S4100-4712
10 holes R	S4100-4713
11 holes R	S4100-4714



Reconstruction Distal Tibial Lateral
Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes L	S4100-4801
6 holes L	S4100-4802
7 holes L	S4100-4803
8 holes L	S4100-4804
9 holes L	S4100-4805
10 holes L	S4100-4806
11 holes L	S4100-4807
5 holes R	S4100-4808
6 holes R	S4100-4809
7 holes R	S4100-4810
8 holes R	S4100-4811
9 holes R	S4100-4812
10 holes R	S4100-4813
11 holes R	S4100-4814




Patella Claw



S	S4100-4901
M	S4100-4902
L	S4100-4903

LC-DCP Femoral Plate
(Use 4.5mm Cortical Screw)



5 holes	S4100-5001
6 holes	S4100-5002
7 holes	S4100-5003
8 holes	S4100-5004
9 holes	S4100-5005
10 holes	S4100-5006
12 holes	S4100-5007
14 holes	S4100-5008

DCP Femoral Plate Broad
(Use 4.5mm Cortical Screw)



5 holes	S4100-5101
6 holes	S4100-5102
7 holes	S4100-5103
8 holes	S4100-5104
9 holes	S4100-5105
10 holes	S4100-5106
12 holes	S4100-5107
14 holes	S4100-5108



LC-DCP Femoral Plate Broad
(Use 4.5mm Cortical Screw)



5 holes	S4100-5201
6 holes	S4100-5202
7 holes	S4100-5203
8 holes	S4100-5204
9 holes	S4100-5205
10 holes	S4100-5206
12 holes	S4100-5207
14 holes	S4100-5208

Proximal Femur Condylus Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes	S4100-5301
6 holes	S4100-5302
7 holes	S4100-5303
8 holes	S4100-5304
9 holes	S4100-5305
10 holes	S4100-5306
12 holes	S4100-5307
14 holes	S4100-5308

Proximal Femur Condylus Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes	S4100-5401
6 holes	S4100-5402
7 holes	S4100-5403
8 holes	S4100-5404
9 holes	S4100-5405
10 holes	S4100-5406
12 holes	S4100-5407
14 holes	S4100-5408



Distal Femoral Medial Plate
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes L	S4100-5501
7 holes L	S4100-5502
9 holes L	S4100-5503
11 holes L	S4100-5504
5 holes R	S4100-5505
7 holes R	S4100-5506
9 holes R	S4100-5507
11 holes R	S4100-5508

Femoral Condylus Buttress Plate-I
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes L	S4100-5601
6 holes L	S4100-5602
7 holes L	S4100-5603
8 holes L	S4100-5604
9 holes L	S4100-5605
10 holes L	S4100-5606
11 holes L	S4100-5607
5 holes R	S4100-5608
6 holes R	S4100-5609
7 holes R	S4100-5610
8 holes R	S4100-5611
9 holes R	S4100-5612
10 holes R	S4100-5613
11 holes R	S4100-5614



Femoral Condylus Buttress Plate-II
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



Femoral Condylus Buttress Plate-III
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes L	S4100-5701
6 holes L	S4100-5702
7 holes L	S4100-5703
8 holes L	S4100-5704
9 holes L	S4100-5705
10 holes L	S4100-5706
11 holes L	S4100-5707
5 holes R	S4100-5708
6 holes R	S4100-5709
7 holes R	S4100-5710
8 holes R	S4100-5711
9 holes R	S4100-5712
10 holes R	S4100-5713
11 holes R	S4100-5714
5 holes L	S4100-5801
6 holes L	S4100-5802
7 holes L	S4100-5803
8 holes L	S4100-5804
9 holes L	S4100-5805
10 holes L	S4100-5806
11 holes L	S4100-5807
5 holes R	S4100-5808
6 holes R	S4100-5809
7 holes R	S4100-5810
8 holes R	S4100-5811
9 holes R	S4100-5812
10 holes R	S4100-5813
11 holes R	S4100-5814



Femoral Condylus Buttress Plate-IV
(Use 4.5mm Cortical Screw,6.5mm
Cancellous Screw)



5 holes L	S4100-5901
6 holes L	S4100-5902
7 holes L	S4100-5903
8 holes L	S4100-5904
9 holes L	S4100-5905
10 holes L	S4100-5906
11 holes L	S4100-5907
5 holes R	S4100-5908
6 holes R	S4100-5909
7 holes R	S4100-5910
8 holes R	S4100-5911
9 holes R	S4100-5912
10 holes R	S4100-5913
11 holes R	S4100-5914



135° DHS Plate
(Use Lag Screw, 4.5mm Cortical
Screw)



3 holes	S4100-6201
4 holes	S4100-6202
5 holes	S4100-6203
6 holes	S4100-6204
7 holes	S4100-6205
8 holes	S4100-6206
9 holes	S4100-6207
10 holes	S4100-6208
11 holes	S4100-6209
12 holes	S4100-6210

Lag Screw



50mm	S4100-6301
55mm	S4100-6302
60mm	S4100-6303
65mm	S4100-6304
70mm	S4100-6305
75mm	S4100-6306
80mm	S4100-6307
85mm	S4100-6308
90mm	S4100-6309
95mm	S4100-6310
100mm	S4100-6311



Reconstruction Plate
(Use 3.5mm Cortical Screw,4.0mm
Cancellous Screw)



4 holes	S4100-6701
6 holes	S4100-6702
8 holes	S4100-6703
10 holes	S4100-6704
12 holes	S4100-6705

Calcaneus Plate-IV
(Use 4.0mm Cancellous Screw)



Long 3 holes	S4100-6601
Long 4 holes	S4100-6602
Short 3 holes	S4100-6603
Short 4 holes	S4100-6604

Calcaneus Plate-III
(Use 4.0mm Cancellous Screw)



Small 10 holes	S4100-6505
Small 12 holes	S4100-6506

Calcaneus Plate-II
(Use 4.0mm Cancellous Screw)



Long 3 holes	S4100-6501
Long 4 holes	S4100-6502
Short 3 holes	S4100-6503
Short 4 holes	S4100-6504



Calcaneus Plate
(Use 4.0mm Cancellous Screw)



S

S4100-6401

M

S4100-6402

Pelvic Plate
(Use 4.0mm Cancellous Screw)



4 holes

S4100-6801

6 holes

S4100-6802

8 holes

S4100-6803

10 holes

S4100-6804

12 holes

S4100-6805

Posterior Wall Plate
(Use 4.0mm Cancellous Screw)



L

S4100-6901

R

S4100-6902

Posterior Column Plate
(Use 4.0mm Cancellous Screw)



L

S4100-7001

R

S4100-7002



Curved Anterior Column Plate
(Use 4.0mm Cancellous Screw)



L

S4100-7101

R

S4100-7102

Angled Anterior Column Plate
(Use 4.0mm Cancellous Screw)



L

S4100-7201

R

S4100-7202

Finger Straight Plate
(Use 2.0mm Cortical Screw)



2 holes

S4100-7301

3 holes

S4100-7302

4 holes

S4100-7303

5 holes

S4100-7304

6 holes

S4100-7305

Finger T-shaped Plate
(Use 2.0mm Cortical Screw)



4 holes

S4100-7401

5 holes

S4100-7402

6 holes

S4100-7403



Finger L-shaped Plate
(Use 2.0mm Cortical Screw)



4 holes L	S4100-7501
5 holes L	S4100-7502
6 holes L	S4100-7503
4 holes R	S4100-7504
5 holes R	S4100-7505
6 holes R	S4100-7506

Finger Curved L-shaped Plate
(Use 2.0mm Cortical Screw)



4 holes L	S4100-7601
5 holes L	S4100-7602
4 holes R	S4100-7603
5 holes R	S4100-7604

Palm Straight Plate
(Use 2.0mm Cortical Screw)



2 holes	S4100-7701
3 holes	S4100-7702
4 holes	S4100-7703
5 holes	S4100-7704
6 holes	S4100-7705



Palm T-shaped Plate
(Use 2.0mm Cortical Screw)



4 holes	S4100-7801
---------	------------

5 holes	S4100-7802
---------	------------

6 holes	S4100-7803
---------	------------

Palm L-shaped Plate
(Use 2.0mm Cortical Screw)



4 holes L	S4100-7901
-----------	------------

5 holes L	S4100-7902
-----------	------------

6 holes L	S4100-7903
-----------	------------

4 holes R	S4100-7904
-----------	------------

5 holes R	S4100-7905
-----------	------------

6 holes R	S4100-7906
-----------	------------

Palm Curved L-shaped Plate
(Use 2.0mm Cortical Screw)



4 holes L	S4100-8001
-----------	------------

5 holes L	S4100-8002
-----------	------------

4 holes R	S4100-8003
-----------	------------

5 holes R	S4100-8004
-----------	------------



2.0mm Cortical Screw(Square Head)
Full-threaded



2.0*6mm	S4100-8101
2.0*8mm	S4100-8102
2.0*10mm	S4100-8103
2.0*12mm	S4100-8104
2.0*14mm	S4100-8105
2.0*16mm	S4100-8106
2.0*18mm	S4100-8107
2.0*20mm	S4100-8108

3.5mm Cortical Screw
Full-threaded



3.5*12mm	S4100-8301
3.5*14mm	S4100-8302
3.5*16mm	S4100-8303
3.5*18mm	S4100-8304
3.5*20mm	S4100-8305
3.5*22mm	S4100-8306
3.5*24mm	S4100-8307
3.5*26mm	S4100-8308
3.5*28mm	S4100-8309
3.5*30mm	S4100-8310
3.5*32mm	S4100-8311
3.5*34mm	S4100-8312
3.5*36mm	S4100-8313
3.5*38mm	S4100-8314
3.5*40mm	S4100-8315
3.5*42mm	S4100-8316



4.5mm Cortical Screw
Full-threaded



4.5*20mm	S4100-8501
4.5*22mm	S4100-8502
4.5*24mm	S4100-8503
4.5*26mm	S4100-8504
4.5*28mm	S4100-8505
4.5*30mm	S4100-8506
4.5*32mm	S4100-8507
4.5*34mm	S4100-8508
4.5*36mm	S4100-8509
4.5*38mm	S4100-8510
4.5*40mm	S4100-8511
4.5*42mm	S4100-8512
4.5*44mm	S4100-8513
4.5*46mm	S4100-8514
4.5*48mm	S4100-8515
4.5*50mm	S4100-8516
4.5*52mm	S4100-8517
4.5*54mm	S4100-8518
4.5*56mm	S4100-8519
4.5*58mm	S4100-8520
4.5*60mm	S4100-8521



4.0mm Cancellous Screw
Full-threaded



4.0*12mm	S4100-8916
4.0*14mm	S4100-8917
4.0*16mm	S4100-8918
4.0*18mm	S4100-8919
4.0*20mm	S4100-8920
4.0*22mm	S4100-8901
4.0*24mm	S4100-8902
4.0*26mm	S4100-8903
4.0*28mm	S4100-8904
4.0*30mm	S4100-8905
4.0*32mm	S4100-8906
4.0*34mm	S4100-8907
4.0*36mm	S4100-8908
4.0*38mm	S4100-8909
4.0*40mm	S4100-8910
4.0*42mm	S4100-8911
4.0*44mm	S4100-8912
4.0*46mm	S4100-8913
4.0*48mm	S4100-8914
4.0*50mm	S4100-8915



4.0mm Cancellous Screw
Half-threaded



4.0*20mm	S4100-8921
4.0*22mm	S4100-8922
4.0*24mm	S4100-8923
4.0*26mm	S4100-8924
4.0*28mm	S4100-8925
4.0*30mm	S4100-8926
4.0*32mm	S4100-8927
4.0*34mm	S4100-8928
4.0*36mm	S4100-8929
4.0*38mm	S4100-8930
4.0*40mm	S4100-8931
4.0*42mm	S4100-8932
4.0*44mm	S4100-8933
4.0*46mm	S4100-8934
4.0*48mm	S4100-8935
4.0*50mm	S4100-8936



6.5mm Cancellous Screw
Full-threaded



6.5*40mm	S4100-9301
6.5*45mm	S4100-9302
6.5*50mm	S4100-9303
6.5*55mm	S4100-9304
6.5*60mm	S4100-9305
6.5*65mm	S4100-9306
6.5*70mm	S4100-9307
6.5*75mm	S4100-9308

6.0mm Double Ends Cannulated
Compression Screw



6.0*60mm	S4100-9801
6.0*65mm	S4100-9802
6.0*70mm	S4100-9803
6.0*75mm	S4100-9804
6.0*80mm	S4100-9805
6.0*85mm	S4100-9806
6.0*90mm	S4100-9807
6.0*95mm	S4100-9808
6.0*100mm	S4100-9809
6.0*105mm	S4100-9810
6.0*110mm	S4100-9811
6.0*115mm	S4100-9812
6.0*120mm	S4100-9813



Herbert Cannulated Compression
Screw



3.0*14mm	/	/	T4100-9601
3.0*16mm	/	/	T4100-9602
3.0*18mm	/	/	T4100-9603
3.0*20mm	/	/	T4100-9604
3.0*22mm	/	/	T4100-9605
3.0*24mm	/	/	T4100-9606
3.0*26mm	/	/	T4100-9607
3.0*28mm	/	/	T4100-9608
3.0*30mm	/	/	T4100-9609
4.0*26mm	/	/	T4100-9610
4.0*28mm	/	/	T4100-9611
4.0*30mm	/	/	T4100-9612
4.0*32mm	/	/	T4100-9613
4.0*34mm	/	/	T4100-9614
4.0*36mm	/	/	T4100-9615
4.0*38mm	/	/	T4100-9616
4.0*40mm	/	/	T4100-9617
4.0*42mm	/	/	T4100-9618
4.0*44mm	/	/	T4100-9619
4.0*46mm	/	/	T4100-9620
4.0*48mm	/	/	T4100-9621
4.0*50mm	/	/	T4100-9622



Cannulated Headless Compression Screw



3.0*12mm	/	/	T4100-9501
3.0*14mm	/	/	T4100-9502
3.0*16mm	/	/	T4100-9503
3.0*18mm	/	/	T4100-9504
3.0*20mm	/	/	T4100-9505
3.0*22mm	/	/	T4100-9506
3.0*24mm	/	/	T4100-9507
3.0*26mm	/	/	T4100-9508
3.0*28mm	/	/	T4100-9509
3.0*30mm	/	/	T4100-9510
3.5*14mm	/	/	T4100-9511
3.5*16mm	/	/	T4100-9512
3.5*18mm	/	/	T4100-9513
3.5*20mm	/	/	T4100-9514
3.5*22mm	/	/	T4100-9515
3.5*24mm	/	/	T4100-9516
3.5*26mm	/	/	T4100-9517
3.5*28mm	/	/	T4100-9518
3.5*30mm	/	/	T4100-9519
3.5*32mm	/	/	T4100-9520
3.5*34mm	/	/	T4100-9521

3.5*36mm	/	/	T4100-9522
3.5*38mm	/	/	T4100-9523
3.5*40mm	/	/	T4100-9524
4.0*16mm	/	/	T4100-9525
4.0*18mm	/	/	T4100-9526
4.0*20mm	/	/	T4100-9527
4.0*22mm	/	/	T4100-9528
4.0*24mm	/	/	T4100-9529
4.0*26mm	/	/	T4100-9530
4.0*28mm	/	/	T4100-9531
4.0*30mm	/	/	T4100-9532
4.0*32mm	/	/	T4100-9533
4.0*34mm	/	/	T4100-9534
4.0*36mm	/	/	T4100-9535
4.0*38mm	/	/	T4100-9536
4.0*40mm	/	/	T4100-9537
4.0*42mm	/	/	T4100-9538
4.0*44mm	/	/	T4100-9539
4.0*46mm	/	/	T4100-9540
4.0*48mm	/	/	T4100-9541
4.0*50mm	/	/	T4100-9542



Cannulated Screw Partially-threaded



3.5*20mm	S4100-9701	\$13.00	T4100-9701
3.5*22mm	S4100-9702	\$13.00	T4100-9702
3.5*24mm	S4100-9703	\$13.00	T4100-9703
3.5*26mm	S4100-9704	\$13.00	T4100-9704
3.5*28mm	S4100-9705	\$13.00	T4100-9705
3.5*30mm	S4100-9706	\$13.00	T4100-9706
3.5*32mm	S4100-9707	\$13.00	T4100-9707
3.5*34mm	S4100-9708	\$13.00	T4100-9708
3.5*36mm	S4100-9709	\$13.00	T4100-9709
3.5*38mm	S4100-9710	\$13.00	T4100-9710
3.5*40mm	S4100-9711	\$13.00	T4100-9711
3.5*42mm	S4100-9712	\$13.00	T4100-9712
3.5*44mm	S4100-9713	\$13.00	T4100-9713
3.5*46mm	S4100-9714	\$13.00	T4100-9714
3.5*48mm	S4100-9715	\$13.00	T4100-9715
3.5*50mm	S4100-9716	\$13.00	T4100-9716
3.5*52mm	S4100-9717	\$13.00	T4100-9717
3.5*54mm	S4100-9718	\$13.00	T4100-9718
3.5*56mm	S4100-9719	\$13.00	T4100-9719
3.5*58mm	S4100-9720	\$13.00	T4100-9720
3.5*60mm	S4100-9721	\$13.00	T4100-9721
4.0*30mm	S4100-9722	\$13.00	T4100-9722

4.0*32mm	S4100-9723	\$13.00	T4100-9723
4.0*34mm	S4100-9724	\$13.00	T4100-9724
4.0*36mm	S4100-9725	\$13.00	T4100-9725
4.0*38mm	S4100-9726	\$13.00	T4100-9726
4.0*40mm	S4100-9727	\$13.00	T4100-9727
4.0*42mm	S4100-9728	\$13.00	T4100-9728
4.0*44mm	S4100-9729	\$13.00	T4100-9729
4.0*46mm	S4100-9730	\$13.00	T4100-9730
4.0*48mm	S4100-9731	\$13.00	T4100-9731
4.0*50mm	S4100-9732	\$13.00	T4100-9732
4.0*52mm	S4100-9733	\$13.00	T4100-9733
4.0*54mm	S4100-9734	\$13.00	T4100-9734
4.0*56mm	S4100-9735	\$13.00	T4100-9735
4.0*58mm	S4100-9736	\$13.00	T4100-9736
4.0*60mm	S4100-9737	\$13.00	T4100-9737
4.5*30mm	S4100-9738	\$13.00	T4100-9738
4.5*32mm	S4100-9739	\$13.00	T4100-9739
4.5*34mm	S4100-9740	\$13.00	T4100-9740
4.5*36mm	S4100-9741	\$13.00	T4100-9741
4.5*38mm	S4100-9742	\$13.00	T4100-9742
4.5*40mm	S4100-9743	\$13.00	T4100-9743
4.5*42mm	S4100-9744	\$13.00	T4100-9744
4.5*44mm	S4100-9745	\$13.00	T4100-9745
4.5*46mm	S4100-9746	\$13.00	T4100-9746
4.5*48mm	S4100-9747	\$13.00	T4100-9747
4.5*50mm	S4100-9748	\$13.00	T4100-9748
4.5*52mm	S4100-9749	\$13.00	T4100-9749
4.5*54mm	S4100-9750	\$13.00	T4100-9750
4.5*56mm	S4100-9751	\$13.00	T4100-9751
4.5*58mm	S4100-9752	\$13.00	T4100-9752
4.5*60mm	S4100-9753	\$13.00	T4100-9753



7.3mm Cannulated Screw
Partially-threaded



7.3*70mm	S4100-9901
7.3*75mm	S4100-9902
7.3*80mm	S4100-9903
7.3*85mm	S4100-9904
7.3*90mm	S4100-9905
7.3*95mm	S4100-9906
7.3*100mm	S4100-9907
7.3*105mm	S4100-9908
7.3*110mm	S4100-9909
7.3*115mm	S4100-9910



Suture Wire



0.4*50m	S7100-9701
0.6*44m	S7100-9702
0.8*12m	S7100-9703
1.0*7.5m	S7100-9704
1.2*5.5m	S7100-9705

Titanium Wire (Ti6Al4V Cable grip)



1.0*520mm, use for patella fracture,tibial fracture and humeral facture			7100-9801
1.3*520mm, use for patella fracture,tibial fracture and humeral facture			7100-9802
1.7*520mm use for femoral fracture and hip protheseses.			7100-9803







INNOVATIVE
ORTHO MEDICAL

1.5mm
Plate